



# DRUG & ALCOHOL TESTING PRODUCTS

***Self-contained Drug Screen Cups: Just add urine and read results in 5 minutes... it's that easy.***

***Features a NEW HINGED LID with key!***

- Accurate • Self-contained, easy to use, no reagents needed
- Stores at room temperature • Rapid results in 5 minutes or less
- Leak-proof cup may be sent to lab for confirmation • Results can be photocopied

## ACCUTEST SPLIT CUPS

*Separate testing area insures integrity of sample*

### ACCUTEST SPLIT CUP 5

- Amphetamine, Cocaine, Opiates, Methamphetamine, THC

Catalog No. DS02AC625

- Amphetamine, Cocaine, Opiates, PCP, THC

Catalog No. DS03AC625

### ACCUTEST SPLIT CUP 6

***with Adulteration Screen***

(pH, Oxidants and Specific Gravity)

- Amphetamine, Cocaine, Opiates, Methamphetamine, PCP, THC

Catalog No. DS08AC625



### ACCUTEST SPLIT CUP 8

***with Adulteration Screen***

(pH, Oxidants and Specific Gravity, Nitrite, Glutaraldehyde and Creatinine)

- Amphetamine, Cocaine, Methamphetamine, Opiates, PCP, THC
- Benzodiazepines, Barbiturates

Catalog No. DS208

### ACCUTEST SPLIT CUP 10

- Amphetamine, Barbiturates, Benzodiazepines, Cocaine, Ecstasy (MDMA), Methadone, Methamphetamine, Opiates PCP, THC

Catalog No. DS07AC625

### ACCUTEST SPLIT CUP 12

- Amphetamine, Cocaine, Marijuana, Methamphetamine, Opiates, Phencyclidine, Benzodiazepines, Methadone, Barbiturates, Oxycodone, Propoxyphene, Ecstasy (MDMA)

Catalog No. DS012





## ACCUTEST DRUG TEST CUP

### 5-PANEL

- Amphetamine, Cocaine, Opiates, PCP, THC [Catalog DS86S625](#)
- Amphetamine, Cocaine, Methamphetamine, Opiates, THC with Adulteration Screen (pH, Oxidants and Specific Gravity) [Catalog DS84](#)

### 6-PANEL

- Amphetamine, Cocaine, Methamphetamine, Opiates, PCP, THC [Catalog DS87S625](#)

### 12-PANEL

- Amphetamine, Barbiturates, Benzodiazepines, Cocaine, Ecstasy, Marijuana, Methadone, Methamphetamine, Opiates, Oxycodone Phencyclidine, Tricyclic Antidepressants [Catalog DS822](#)



## ACCUTEST MULTIDRUG SCREEN

- Amphetamine, Cocaine, Opiates, Methamphetamine, THC [Catalog No. DS01AC425](#)
- Amphetamine, Cocaine, Opiates, PCP, THC [Catalog No. DS61AC425](#)
- Rapid results in 5 minutes or less
- Available in single or multiple combinations
- CLIA waived cassette format also available



## INTOXILYZERS

**AUTOMATIC SAMPLING** handheld breath alcohol screeners

### INTOXILYZER 400

- DOT Screening
- Results displayed & printed
- Fully automatic
- Designed for employee testing programs

[Catalog No. DS400](#)

### INTOXILYZER S-D5

- Last test recall
- Fuel cell accuracy
- Large LED display
- Easy to use

[Catalog No. DS45CM501](#)



## ACCUSTRIP SALIVA ALCOHOL TEST

- One-step rapid saliva test quickly detects the presence of alcohol
- Provides results in 2 minutes • Strips Individually Packaged - 25 tests/box
- Detection range is from 0.02% to 0.30% for approximate relative blood alcohol concentration • Accurate, reliable and easy to use
- Includes convenient collection cup
- Room temperature storage • **FDA-cleared, CLIA-waived** [Catalog No. DS97AC125](#)



## ACCUTEST BREATH ALCOHOL DETECTOR

- One Step, rapid, semi-quantitative breath alcohol test
- Indicates relative Blood Alcohol Concentration (BAC) at 0.02% cutoff
- Results in 2 minutes • Individually Packaged - 20 tests/box
- Accurate and Reliable
- Easy to Use • Room temperature Storage
- **FDA-cleared for OTC Use** [Catalog No. DS95AC220](#)



### JANT PHARMACAL CORPORATION

16530 Ventura Boulevard, Suite 512, Encino CA 91436  
800.676.5565 818.986.8530 Fax 818.986.0235  
[www.jantdx.com](http://www.jantdx.com) [info@jantdx.com](mailto:info@jantdx.com)

Accutest is a registered trademark of Jant Pharmacal Corporation

