



H. PYLORI UREASE TEST CASSETTE



For Rapid Detection of Helicobacter Pylori in Gastric Biopsies

ACCUTEST H. PYLORI UREASE TEST FEATURES

NO REFRIGERATION REQUIRED	Can be stored at room temperature.
SIMPLE & CONVENIENT	May be used in any CLIA-waived endoscopy unit with minimal training. Simply place biopsy in the test well, close and wait for results.
RAPID RESULTS	Positive results in as little as 1 minute, negative results in 1 hour.
ACCURATE	>99% Sensitivity and Specificity.*
DRY TEST	No incubation or reagents required.
COST EFFECTIVE	Low cost per test.
3 YEAR SHELF LIFE	Good stability, superior shelf life.
CPT CODE	87077QW



*Refer to the Instructions for Use (IFU) for detailed information and studies.


CATALOG #	DESCRIPTION	UNITS
ID750	Accutest H. pylori Urease Test	50 tests/box (control included)
IDC73	Accutest Urease/Glycerin Controls Solution	3 mL dropper bottle

JANT PHARMACAL CORPORATION
LABORATORY & POINT OF CARE DIAGNOSTICS

accutest[®] is a registered trademark of Jant Pharmacal Corporation.

 800.676.5565

 info@jantdx.com

 jantdx.com